



0.1µm AseptiCap KS

Sterilization Grade

Hydrophilic Polyethersulfone (PES) Membrane Devices for Liquid Streams in Biopharmaceuticals

Data Sheet

Biopharmaceutical processing requires sterilizing grade microfiltration at multiple stages to meet specific process requirements.

Processes managers are continuously looking for microfiltration solutions to upstream, downstream, intermediate processes and final biological preparations. Since bio manufacturing is a multi stage process and bio molecules by nature are extremely sensitive, they are looking for:

- Minimizing protein losses due to adsorption to improve over all product yields
- Minimizing filter extracts which add up due to multiple points of use in a process
- High throughputs to achieve process economy
- Choice of filter end connections for easy and reliable quick connections
- > Absolute retentions for higher sterility assurance

mdi produces a wide range of Sterilizing grade PES membrane devices to meet filtration requirements of biopharmaceutical processing. These filter devices are validated to meet compendia and regulatory requirements and are well characterized. They meet key process requirements such as high retention efficiency, very high protein recoveries, extremely low extractables, high throughputs, wide chemical compatibility and other important characteristics.

With the advantages of pre filtration layer built into the device for higher throughputs, linear scalability of filter area for smooth transitions from lab scale to pilot to process scale and widest range of end connections for quick and reliable connections to the existing fittings, **mdi** *AseptiCap* KS filters are a universal solution for process filtration.

AseptiCap KS

Datasheet

PES Membrane Devices

for Biopharmaceuticals

Asepticap KS 0.1 micron capsule filters uses **mdi** PES membrane in Polypropylene housing. No adhesives or glue are used in the manufacturing process and all bonding is done by heat welding.

The products are deeply validated for use in Biopharmaceutical applications and specially recommended for single use systems. *Asepticap KS* are manufactured in class 10,000 clean rooms and ISO 9001 certified facilities.

Applications

Sterile Filtration of

- Cell culture media
- > Cell culture media containing serum
- > Media additives
- pH adjusters
- > Final product concentrates

Key Features

- Absolute retention
- > 100% integrity tested
- Low protein binding
- > Very low hold up volume in filters
- > High flow rates
- Serial construction with prefilter for higher throughput with fouling streams
- Bioburden maintained below 1000 cfu/device
- Endotoxin level certified to be <0.25 EU/ml</p>
- Widest range of end connections
- Products available for total scalability from a few ml to thousands of liters
- Total traceability through unique serial number for each filter
- > Individual certificate of quality for each device
- Sterilizable by EO gas or autoclaving

Validation Services

The regulatory requirements emphasize on the need to validate the efficacy of the 'Sterilizing Filter' with drug product under simulated worst-case conditions of use.

mdi provides validation services supported by customized validation protocols and world class test facilities to assist you in filter validations with your specific drug product.

Quality Assurance

mdi's quality management system emphasizes on quality by design rather by end product testing. Robust processes are developed for product manufacturing and are continuously monitored to ensure that the products meet their predetermined specifications and lot to lot reproducibility is ensured.

Certificate of Quality

Each capsule filter is accompanied by individual certificate of quality to ensure traceable documentation at user's end.

It certifies the product compliance to various regulatory as well as user requirements.

Validated for Microbial Retention

Integrity test data have been correlated to actual microbial retention with *Acholeplasma laidlawii* ATCC 23206 at a challenge level $\geq 10^7$ organisms per cm² to establish acceptable integrity test values. Also validated for retention of *B. diminuta* ATCC 19146 as per ASTM F838-05.

100% Integrity Tested

Each *AseptiCap KS* is tested for integrity to comply with validated Acceptable Integrity Test Specifications.

Flow Rate

Each lot is tested for clean water flow rates to ensure that flow rates are within the specifications.

Adsorption

AseptiCap KS filters are validated for low protein binding to ensure minimal active ingredient losses when used for filtration of high value proteins.

Pressure, Temperature Endurance

AseptiCap KS filters are validated to endure high operating pressure and temperature conditions which may be encountered during use.

These filters are also validated for high burst pressure to ensure user safety in case of inadvertent pressure build-up.

Extractables

Extractables/leachables from sterilizing filters, used at various stages of a biopharmaceutical manufacturing process, will add on and may impact the impurity profile of the desired product.

AseptiCap KS filters are validated to exhibit low extractables under harsh extraction conditions.

Bioburden Testing

Device bioburden is tested as per ISO 117 37-1 and assured to be <1000 cfu/device.

Endotoxin Testing

Aqueous extracts exhibit < 0.25 EU/ml as established by Limulus Amebocyte Lysate (LAL) Test

Total Traceability

AseptiCap KS filters come with completely traceable lot numbers and unique identification number to facilitate easy and fast retrieval of manufacturing and quality control data.

These unique lot and identification numbers are laser etched on each filter device and also printed on the labels of the box in which individual filter is packed.

Packaging Integrity

AseptiCap KS filters are fitted with vent caps and are packed in pouch to ensure package integrity during transit as well as to prevent particulate contamination while transferring to clean room assembly or process areas.

Other Regulatory Compliance

- Complies with USFDA 21 CFR 210.3(b)(6) for fiber release
- Complies with USFDA 21 CFR 177.1520 for fractional dissolution
- Materials of construction tested for toxicity as per Biological Reactivity Tests, In vivo, USP <88> for class VI Plastics
- Complete filter devices tested for cytotoxicity as per Biological Reactivity Tests, In vitro, USP <87>

Performance Data


Datasheet

Low Protein Binding

A comparative study on **mdi** PES membrane exhibits much lower protein adsorption than other competing membranes of Cellulose Acetate and PVDF.


The low protein binding results in increased overall product yield and higher throughputs with biological streams.

Protein Binding (µg/cm²)



0.1 μm <i>AseptiCap</i> Filters	Protein Binding
25 mm, 5 cm ²	1.7 µg
50 mm, 20 cm ²	7 µg
1″, 250 cm²	88 µg
2", 500 cm ²	187 µg
10″, 6000 cm²	2275 µg

High Throughputs



The high throughput translates to lower filtration costs, less number of filter changes and overall economy of operations.

Very Low Hold-Up Volumes

mdi PES membrane capsule filters are designed to offer very low hold-up volumes to minimize filtration losses and maximize product recovery.

Filter Devices	EFA* (Nominal)	Hold up Volume
AseptiCap KS, 25mm	5cm²	< 50µl
AseptiCap KS, 50mm	20cm ²	< 200µl
AseptiCap KS, 1"	250cm ²	< 5ml
AseptiCap KS, 2"	500cm ²	< 25ml
AseptiCap KS, 5"	1000cm ²	< 45ml
AseptiCap KS, 8"	2000cm ²	< 60ml

*EFA: Effective Filtration Area

Performance Data

Datasheet

Extractables


It is useful to evaluate extractables that may be leeched out of the filter and enter the process stream. **mdi** filters give low extractables under harsh preconditioning and extraction conditions.

Low extractables mean less addition to impurity profile of the biological product from the filters.


Extraction Time: 24 hours

	Non Volatile Residue							
Model Solvent	AseptiCap KS 1" (250 cm²)	AseptiCap KS 10" (6000 cm²)						
Water @ 35 °C	1.6 mg	38.26 mg						
Water @ 80 °C	1.8 mg	43.04 mg						

FTIR Analysis of Extractables From *AseptiCap KS* 1" Capsule Filter with Water @ 35 °C




FTIR Analysis of Extractables From *AseptiCap KS* 1" Capsule Filter with Water @ 80 °C




	Non Volatile Residue						
Model Solvent	<i>AseptiCap K</i> S 1″ (250 cm²)	AseptiCap KS 10" (6000 cm²)					
100% Ethanol @ 35 °C	13.4 mg	320.43 mg					

FTIR Analysis of Extractables From AseptiCap KS 1" Capsule Filter with 100% Ethanol @ 35 °C



FTIR Analysis of Extractables From *AseptiCap KS* 1" Capsule Filter with 100% Ethanol @ 80 °C



The Spectrum of extracts from *AseptiCap KS* capsule filters with 100% ethanol under extreme extraction conditions show presence of various components used in the manufacture of **mdi** PES membrane capsule filters.

Easy Connect

Datasheet

Widest Range of End Connections

Biopharmaceutical processes involve transfer of high value fluids through multiple process steps. Making high quality, reliable, flexible and functionally convenient connectivity with filters is of utmost value to the bio-processors.

mdi AseptiCap KS filters offer a wide range of reliable end connections for functional convenience and customized connectivity.

Validated for Performance

These end connections are manufactured with tight dimension tolerance and are validated for strength and connection integrity under extreme use conditions as well as for their ability to withstand prevalent sterilization methods including EO sterilization and autoclaving.

Customized Connectivity

mdi AseptiCap KS filters are available in a wide range of end connections and are also customized to offer different inletoutlet combinations to meet the unique connectivity needs in biopharmaceutical process assemblies where, for example, stainless steel components with sanitary flange connections are sometimes required to be connected to single use disposable systems through quick-connectors or hose barb connections.



1¹/₂" Sanitary Flange to ¹/₂"Barb Hose



³⁄₄" Sanitary Flange



½″ HB



1⁄4″ SHB



1¹/₂" Sanitary Flange



1/2" Single Stepped HB



Quick Connector

Some end connections available with AseptiCap

1½" Sanitary Flange to ¾" Sanitary Flange





AseptiCap with HighSecurity ¹/₂" hose barb connection

Linear Upscaling from R&D to Production Process

Datasheet

Scientists are concerned about filter fluid interaction impacting the stability, purity, strength etc. of the drug product, and they take a keen interest in filter selection at the formulation development stage itself. Although preliminary compatibility data support initial filter selection, for stability studies detailed filter validations are required to provide enough documented evidence to justify specific filter use.

A critical requirement that needs to be addressed at this stage is of scalability from R&D to pilot scale to full scale production processes.

mdi offers a wide range of *AseptiCap KS* filters to provide linear scale up from lab scale to production process. While scaling up the process, the appropriate size filter can be selected by increasing the effective filtration area of filter proportionate to the process fluid volumes.

All Materials of construction as well as manufacturing process are identical for all filter devices starting from 5 cm² to 18000cm² hence process scaling can be facilitated without triggering additional validation studies for given process conditions. **mdi** provides complete documentation for each of the *AseptiCap KS* filters there by reducing the additional validation cost and time.



AseptiCap KS 25mm, 5cm²



AseptiCap KS 50mm, 20cm²



AseptiCap KS 1", 250cm²



AseptiCap KS 2", 500cm²



AseptiCap KS 5", 1000cm²



AseptiCap KS 8", 2000cm²

Filter Devices	EFA* (Nominal)	Hold up Volume
AseptiCap KS, 25 mm	5cm ²	< 50µl
AseptiCap KS, 50 mm	20cm ²	< 200µl
AseptiCap KS, 1"	250cm ²	< 5ml
AseptiCap KS, 2"	500cm ²	< 25ml
AseptiCap KS, 5"	1000cm ²	< 45ml
AseptiCap KS, 8"	2000cm ²	< 60ml
AseptiCap KS, 10"	6000cm ²	-
AseptiCap KS, 20"	12000cm ²	-
AseptiCap KS, 30"	18000cm ²	-



AseptiCap KS 10", 6000cm²

Specifications

Datasheet

0.1 µm *AseptiCap KS* (with Prefilter)

		Construction						
Membrane		0.1 μm Hydro	philic PES					
Prefilter Mem	brane	0.2 μm or 0.45 μm Hydrophilic PES						
Plastic parts		Polypropylene						
		Integrity Testing/ Retention						
Bubble Point		\geq 31 psi (2.18 Kg/cm ²) with 50% IPA						
		LRV> 7 for Acholeplasma laidlawii ATCC 23206 p	per cm ² of filter area					
Bacterial Rete	ention	LRV> 7 for Brevudimonas diminuta ATCC 19146	per cm ² of filter area as per ASTM F 838-05					
		Size						
Size		25mm	50mm					
Effective Filtra	ation Area (Nominal)	5 cm ²	20 cm ²					
	1⁄4″ SHB I/O	_	79 mm					
Dimensions	, · · · · · · · · · · · · · · · · · · ·	-	51 mm					
(End to End)	Female Luer Lock Inlet/	23 mm						
Male Luer Slip Outlet Operational Radius (with Vent/ Drain)		15 mm	28 mm					
operationari		Operational	20 1111					
Max Operatiu	ng Temperature	55 ℃	60 °C					
Max. Operating Temperature Max. Differential Pressure		75 psi (5 Kg/cm ²) @ 25 °C	42 psi (3 Kg/cm²) @ 30 °C					
Max. Differen	By Gas	Sterilizable by Ethylene Oxide						
Sterilization	By Autoclave							
	,	Autoclavable at 125 °C for 30 minutes, 25 Cycles. Can not be in-line steam sterilized						
Shelf Life		3 year after EO sterilization						
		Assurance						
Toxicity		Passes Bioreactivity test, In Vivo, as per USP <8	8> for Class VI plastics					
Cytotoxicity		Passes Biological Reactivity Tests, In vitro, USP <87> for cytotoxicity						
Bacterial Endo	otoxin	Aqueous extracts exhibit < 0.25 EU/ml as established by Limulus Amebocyte Lysate (LAL) Test as per USP <85>						
Non Fiber Rele	easing	Passes test as per USP and comply with USFDA 21 CFR Part 210.3(b)(6) for fiber release						
TOC and Cond	ductivity	Meets the WFI requirements of USP for TOC <643> and Conductivity <645> after a specified minimal flush						
pH Compatibi	ility	Compatible with pH range of 1 - 10						
Extractables v	vith WFI	Passes NVR test as per USP <661>						
Indirect Food	Additives	Comply with USFDA 21 CFR Part 177.1520						
Oxidizable Su	bstances	Passes test as per USP <1231>						
Quality Mana	gement System	ISO-9001 Certified						
		DMF No. 015554						

Specifications 0.1µm *AseptiCap KS* (with Prefilter)

Datasheet

		Co	onstruction							
Membrane		0.1 μm Hydrophilic PES								
Upstream Me (in case of As		0.2 μm or 0.45 μm Hydrophilic PES								
Support Laye	ers		Polyest	ter						
Plastic parts			Polyprop	ylene						
		Integrity 1	Testing/ Retention							
Bubble Point	:	<u>></u> 31psi with 50% IPA/	Water solution							
		LRV> 7 for Acholeplasn	na laidlawii ATCC 23206 p	er cm ² of filter area						
Bacterial Ret	ention	LRV> 7 for <i>Brevudimonas diminuta</i> ATCC 19146 per cm ² of filter area as per ASTM F 838-05								
		Size								
Size		1″	2″	5″	8″					
Effective Filt	ration Area (Nominal)	250cm ²	500cm ²	1000cm ²	2000 cm ²					
	Flow Rate @ 10 psi with c Connection	1.3 lpm	2.5 lpm	4.8 lpm	6.5 lpm					
	1½" Sanitary Flange I/O	91 mm	110 mm	161 mm	211 mm					
Dimension	¹ / ₂ " Hose Barb I/O	90 mm	112 mm	164 mm	215 mm					
Dimensions (End to End)	1½" Sanitary Flange Inlet ½" Single Step Hose Barb Outlet	-	111 mm	162 mm	212 mm					
³ 4" Sanitary Flange I/O		91 mm	103 mm	155 mm	205 mm					
Operational F	Radius (with Vent/ Drain)	30 mm	65 mm	65 mm	65 mm					
Vent and Dra	in	1/4" Hose Barb with Silicone "O" rings								
		Operational								
Max. Opera	ting Temperature	80 °C @ < 30 psi (2 Kg/cm²)								
Max. Differe	ential Pressure	60 psi (4 Kg/cm²) @ 30 °C								
Sterilizatior	By Gas	Sterilizable by Ethylene Oxide								
Stermzation	By Autoclave	Autoclavable at 125 °C for 30 minutes, 25 Cycles. Can not be in-line steam sterilized								
Shelf Life		3 year after EO sterilization								
			Assurance							
Toxicity		Passes Bioreactivity te	est, In Vivo, as per USP <88	3> for Class VI plastics						
Cytotoxicity		Passes Biological Reactivity Tests, In vitro, USP <87> for cytotoxicity								
Bacterial End	otoxin	Aqueous extracts exhibit < 0.25 EU/ml as established by Limulus Amebocyte Lysate (LAL) Test as per USP <85>								
Non Fiber Re	leasing	Passes test as per USP	and comply with USFDA	21 CFR Part 210.3(b)(6) fo	r fiber release					
TOC and Con	ductivity	Meets the WFI require	ements of USP for TOC <64	43> and Conductivity <64	5> after a 3 liter flush					
pH Compatik	bility	Compatible with pH r	ange of 1 - 10							
Extractables	with WFI	Passes NVR test as per	r USP <661>							
Indirect Food	l Additives	Comply with USFDA 2	21 CFR Part 177.1520							
Oxidizable Su	ubstances	Passes test as per USF	°<1231>							
	agement System	ISO-9001 Certified								
Quality Mana	igement System	ISO-9001 Certified								

Specifications 0.1µm *AseptiCap KS* (with Prefilter)


Datasheet

		Con	struction						
Membrane			0.1 µm Hydropł	nilic PES					
Upstream Mer (in case of Asep		0.2 μm or 0.45 μm Hydrophilic PES							
Support Layer	S		Polyeste	er					
Plastic parts			Polypropyl	lene					
		Integrity Te	esting/ Retention						
Bubble Point		≥ 31psi with 50% IPA/	Water solution						
		LRV> 7 for Acholeplasn	na laidlawii ATCC 23206 p	per cm ² of filter area					
Bacterial Reter	ntion	LRV> 7 for Brevudimon	as diminuta ATCC 19146	per cm ² of filter area as p	er ASTM F 838-05				
			Size						
Size		5″	10″	20″	30″				
	tion Area (Nominal)	3000cm ²	6000cm ²	12000cm ²	18000 cm ²				
	low Rate @ 10 psi with								
	lange Connection	8 lpm	17 lpm	29 lpm	45 lpm				
Max. Air Diffus (@ 50psi (3.51	Kg/cm ²) with water)	≤ 15 ml/min	≤ 29 ml/min	≤ 58 ml/min	≤ 87 ml/min				
Dimensions	½" Single Step Hose Barb I/O	217 mm	332 mm	607 mm	882 mm				
Dimensions (End to End) Inline Capsule Filters	1½" Sanitary Flange Inlet ½" Single Step Hose Barb Outlet	203 mm	203 mm 332 mm 607 mm		882 mm				
	1½" Sanitary Flange I/O	207 mm	07 mm 326 mm 601 mm						
Operational Ra	dius (with Vent/ Drain)	78 mm 78 mm 78 mm 78							
Vent and Drain	1	1/4" Hose Barb with Sil	icone "O" rings						
		0	perational						
Max. Operati	ng Temperature	80 °C @ < 30 psi (2 Kg/	cm²)						
Max. Differen		60 psi (4 Kg/cm ²) @ 30 °C							
	By Gas	Sterilizable by Ethylene Oxide							
Sterilization	By Autoclave	Autoclavable at 125 °C for 30 minutes, 25 Cycles. Can not be in-line steam sterilized							
Shelf Life		3 year after EO sterilization							
		Α	ssurance						
Toxicity		Passes Bioreactivity te	st, In Vivo, as per USP <88	3> for Class VI plastics					
Cytotoxicity		Passes Biological Reactivity Tests, In vitro, USP <87> for cytotoxicity							
Bacterial Endo	toxin	Aqueous extracts exhibit < 0.25 EU/ml as established by Limulus Amebocyte Lysate (LAL) Test as per USP <85>							
Non Fiber Rele	asing	Passes test as per USP	and comply with USFDA	21 CFR Part 210.3(b)(6) f	or fiber release				
TOC and Cond	luctivity	Meets the WFI require	ments of USP for TOC <64	43> and Conductivity <6	45> after a 3 liter flush				
pH Compatibil	lity	Compatible with pH range of 1 - 10							
Extractables w	ith WFI	Passes NVR test as per USP <661>							
Indirect Food /	Additives	Comply with USFDA 2	1 CFR Part 177.1520						
Oxidizable Sub	ostances	Passes test as per USP	<1231>						
Quality Manag	jement System	ISO-9001 Certified							
USFDA		DMF No. 015554							

Ordering Information

Datasheet

0.1 µm AseptiCap KS 25mm PES Membrane Capsule filter



0.1 µm AseptiCap KS 50mm PES Membrane Capsule filter

Туре		Size		ze Pore Size		Inlet/Outlet		X	х	Sterility		Pack Size	
	Code		Code		Code		Code				Code		Code
AseptiCap KS		50mm	10	0.1µm	36	1⁄4″ SHB	В			Non Sterile	1	12	08
(0.45 µm Upstream)	VKSX					³ ⁄ ₄ " Sanitary	S			EO Sterile	2		
<i>AseptiCap KS</i> (0.2 μm Upstream)	VKS1					Flange			1				
Example:	1		7	,		V		V	V		1	,	V
VKSX		1	0	30	5	BB		x	х	1	1	0	8

Note: Inlet/Outlet Connections and Pack Sizes available with different diameter filters as follows:

Connections Available										
Inlet/Outlet	25mm	50mm								
1/4" - 3/4" Stepped Hose Barb	x	\checkmark								
3/4" Sanitary Flange	x	\checkmark								
Female Luer Lock	Inlet Only	х								
Male Luer Slip	Outlet Only	х								
1/8" Hose Barb		х								
Male Luer Lock	Outlet Only	х								
¼" Hose Barb	\checkmark	х								

Pack Size Available									
Pack Size	25mm	50mm							
12/Pack	х	\checkmark							
100/Pack	\checkmark	х							

Ordering Information

Datasheet

0.1 µm AseptiCap KS PES Membrane Capsule filter

Туре		S	Size Pore S		Size	Inlet/Outlet		Inlet/Outlet		В	ell	Sterilit	y	Pack	Size
	Code		Code		Code			Code			Code		Code		Code
AseptiCap KS		1″	51	0.1µm	36	1/	4″ SHB	A		Yes*	В	Non Sterile	1	1	01
(0.45 µm Upstream)	DKSX	2″	52			½″ H	lose Barb	D		No Bel	Х	EO Sterile	2		
AseptiCap KS		5″	53			1½" San	itary Flange	E		Bell with co	ver C				
(0.2 µm Upstream)	DKS1	8″	57			¾" San	itary Flange	S							
						Quick	Connector	J							
						1/2" Single S	Step Hose Barb	Q	*Bell is available with						
						Female luer lock U			- ½"HB outlet connections in 1", 2", 5" and 8" capsule filters						lters
						Male	e luer slip	W	- ¼" SHB outlet connection in 1" capsule filters only						
						3/16″	Hose Barb	N							
Example:	Example:				3/8″	Hose Barb	I								
DKSX			57			36	DD		X		Х	1		0	1

Note: Inlet/Outlet Connections available with different Sizes/Length as follows:

Inlet/Outlet	Size/Length							
iniet/Outlet	1″	2″	5″	8″				
¼" Stepped Hose Barb	\checkmark		\checkmark	\checkmark				
1/2" Single Step Hose Barb	Х		\checkmark	\checkmark				
1/2"Hose Barb	\checkmark		\checkmark	\checkmark				
1½" Sanitary Flange	\checkmark	\checkmark	\checkmark	\checkmark				
¾" Sanitary Flange	\checkmark	\checkmark	\checkmark	\checkmark				
Quick Connector	\checkmark		\checkmark	\checkmark				
Female Luer Lock	\checkmark	\checkmark	\checkmark	\checkmark				
Male Luer Slip	Outlet Only	х	х	Х				
3/16" Hose Barb	\checkmark	\checkmark	\checkmark	\checkmark				
3/8" Hose Barb	Х	\checkmark	\checkmark	\checkmark				

Bell at Outlet Available with (Size/Outlet)						
1"/ ¼" SHB						
1", 2", 5", 8"/ ½" HB						

0.1 μm AseptiCap KS PES Membrane Large Capsule filter

Туре		Size		Pore Size		Inlet/Outlet		X	Inline/T-Line		Sterility		Pack Size	
	Code		Code		Code		Code			Code		Code		Code
AseptiCap KS		5″	53	0.1µm	36	1/2" Single Step Hose Barb	Q		Inline	Х	Non Sterile	1	1	01
(0.45 µm Upstream)		10″	54			1½" Sanitary Flange	E		T-Line*	Т	EO Sterile	2		
AsentiCan KS	ceptiCap KS LKS1 30"	20″	55			, ,	-							
(0.2 µm Upstream)		30″	56			3/8" Hose Barb	1							
(0.2 µm opstream)						1" Hose Barb	Z							

Example:

LKSX 54	36	EE	Х	т	1	01
---------	----	----	---	---	---	----

*T-line is not available in 5" Capsule filter

*T-line Capsule filter are available with $1\frac{1}{2}$ " Sanitary Flange I/O connection only

Note: Inlet/Outlet Connections available with different Sizes/Length as follows:

Inlet/Outlet		Inli	ine	T-Line			
iniet/Outlet	5″	10″	20″	30″	10″	20″	30″
1/2" Single Step Hose Barb	\checkmark	\checkmark	\checkmark	\checkmark	Х	Х	Х
1 ¹ / ₂ " Sanitary Flange				\checkmark			
3/8" Hose Barb	\checkmark	\checkmark	\checkmark	\checkmark	Х	Х	Х
1" Hose Barb	Х	\checkmark	\checkmark	\checkmark	Х	Х	Х

Advanced Microdevices Pvt. Ltd.

20-21, Industrial Area, Ambala Cantt-133 006, INDIA Tel : +91-171-2699290, 2699471 E-mail : info@mdimembrane.com Website : www.mdimembrane.com

DST DKLK36X1129C