
FilterCheck 06 Integrity Test Equipment



Integrity Testing of Critical Sterilizing Filters



BRC FC06XXX1519H

FilterCheck 06

Integrity testing of sterilizing grade filters in aseptic processes is critical as undetected filter integrity failure may compromise drug product sterility and have huge regulatory and financial implications.

Process engineers are looking for accurate and reliable non destructive integrity testing of sterilizing filters. Online integrity testing of sterilizing filters is even more demanding as integrity test equipment are required to be aseptic area compatible and should be easy to move from one filter location to another. The high incidence of sterilizing filters in a pharmaceutical or biopharmaceutical process requires faster test times in order to meet process demands.

Current testers are unable to meet one or the other need. **mdi** FilterCheck 06 fulfils all the above needs.


FilterCheck 06 is a fully automatic, compact, 21 CFR part 11 compliant integrity tester for reliable testing of your critical filters.

It uses the latest mass flow measurement technology for real time flow measurement to deliver accurate and precise test results. The user friendly, easy to use validated software and calibrated sensors ensure test reliability.

SS 316 L body and a light weight compact design, allows it to be stationed within the aseptic process area for online integrity testing of sterilizing filters. It then wirelessly transfers test results to the laboratory network as well as to a printer.

mdi FilterCheck 06 with fast integrity testing (~15 minutes/test) delivers high productivity. It is supported by extensive documentation including COQ, IQ, OQ, OQ software, PQ, calibration report and user manual.


Unique Performance Advantages



Compatibility with Sterile Processing Areas

- Electropolished SS 316L Surface
- No Instrument Cooling Fan
- Wireless and Paperless to minimize particle generation in sterile areas





Easy Portability

- Smallest Integrity Tester on the Market
- Weight: 9kg



All types of Integrity tests


For Hydrophilic as well as Hydrophobic FiltersOption for QuickTest for Higher Productivity


Main Menu 21 CFR part 11 Compliance Administrator for Data Reliability and Security Supervisor Operator Three Access levels with Password Security System AuditTrail Membrane Technologies System Calibration System Audit Trail Mon, 6 Apr 2015 12:11 PM 4:37:17 PM Current Time : **Self Diagnosis** Time Elapsed : 2:14 \bigcirc ternal Communication nlet Pressure Pressure Calibration Internal Leakages nternal Leakage External Leakages to ensure Test Filter Connectivity Abort Test e · Buhhle Poin arator Name : Admin /in. Bubble Point : 50 psi Shows Real Time Test Status

Shows Near time test Status

Graphical representation of test results

PDF report generation





Inbuilt Equipment Safety with separate Hydrophobic Filters for

- Entrained moisture in Inlet Air
- Inadvertent back flow of test fluids



Wi-Fi Connectivity and Data Back Up

FilterCheck 06 connects with the external printer (with inbuilt Wi-FiRouter) and a remote back up data repository.

The remote back up data repository for test reports can be set up with an optional software specially designed for reliability and security.



Validation and Calibration

FilterCheck 06 uses a completely validated software.

The pressure and mass flow sensors are calibrated with instruments traceable to national reference standards.

Each instrument is supplied with a calibration report.

Equipment Supplied

- > FilterCheck 06: Main Integrity Test Unit
- > Outlet tubing connections with:
 - SS Quick Connector to connect to the main test unit, and 1.5"TC connector to connect to the test filter inlet
 - SS Quick Connector to connect to the main test unit, and Staubli Female (RBE03) connector to connect to the test filter housing
 - SS Quick Connector to connect to the main test unit, and 1/8" Stepped Hose Barb connector to connect to the test capsule filter inlet
- Laser Printer with inbuilt Wi-Fi router
- Power Cord

mdi Membrane Technologies INC.

5340 Jaycee Avenue, Suite A, Harrisburg, PA 17112 Phone: +1-717-412-0943, Fax: +1-717-695-9637 E-mail: rs@mdimembranetech.com Website: www.mdimembranetech.com

Specifications

Dimensions		
Dimensions(Main Test Unit)	35cm X 32cm X 12cm	
Weight	9 Kg	

Electrical data		
Voltage	100-240V AC	
Input Frequency	50 Hz	
Battery Pack	14.4V - 4.5 AH	
Battery Charger (in-built)	Yes	
Power Backup	6 hours	
Operating Range		
Minimum Inlet Pressure	<u>≥</u> 70psi	
Maximum Inlet Pressure	100psi	
Operating Temperature	5°-35° C	
Storage Temperature	8°-30° C	
Relative Humidity	10-80 %	

Measuring Range	
Diffusion Flow	0.2-100 ml/min
Bubble Point	15-80psi
Water Intrusion Test	0.03-10 ml/min

Accuracy			
Diffusion Flow	$\pm 5\%$ or ± 1 ml whichever is greater		
Bubble Point	<u>+</u> 0.7psi		
Water Intrusion Test	$\pm 5\%$ or ± 0.02 ml whichever is greater		

Report Printing

External Laser Printer with Inbuilt Wi-Fi Router

Support Documents

FilterCheck 06 comes with full documentation which is easy to follow.

Certificate of Quality		Yes
Calibration Report		Yes
User Manual		Yes
IQ		Yes
OQ	Test Equipment	Yes
	Application Software	Yes
PQ		Yes

ADVANCED MICRODEVICES PVT. LTD.

20-21, Industrial Area, Ambala Cantt- 133 006, India Tel: +91 - 171-2699290, 2699471 E-mail: info@mdimembrane.com Website: www.mdimembrane.com