

Certificate of Quality

The Polyethersulfone Membrane Capsule filter has been manufactured in a **mdi** facility in compliance with **ISO 9001** regulations using *validated production processes*.

AseptiPrime KS-y PES Membrane Capsule Filter

Catalog No. : DKX15736DDXX301

Type : DK-S

Pore Size : $0.1 \mu m (0.2 \mu m + 0.1 \mu m)$

Lot Number : DK5577G Sl.No. 004

Ster. No. : R002 Expiry Date : 2019 - 07

SPECIFICATION

Length	8"
Filter Media	Polyethersulfone Membrane
Drainage Layers	Polyester
Housing	Polypropylene
Differential Pressure	< 4Kg/cm² at 30 °C
Maximum Operating Temperature	80 °C @ < 2 Kg/cm ²
Sterilization	Pre sterilized by Gamma Irradiation

LOT RELEASE CRITERIA

100% Integrity Tested : The Capsule filter has been tested for integrity by Air Diffusion Flow Test using

purified water and Bubble Point Test using 50% IPA/Water solution.

Diffusion flows was: ≤ 13.5 ml/min @ 3.52 kg/cm²

Bubble point was: \geq 31 psi (2.14 Bar).

Water Flow Rate :≥ 7.5 lpm @ 0.70 Kg/cm² @ 27 °C

Sterility : Samples passed the sterility test in accordance with USP <71>.

VALIDATED FOR

Extractable : Within limits as specified in USP.

Oxidizable matter : Passes test as per USP.

Biosafety : Passes Biological Reactivity Tests, *In Vivo* for Class VI plastic as described

in USP <88>.

Cytotoxicity: Passes Biological Reactivity Tests, *In Vitro* as described in USP <87>.

indirect Food Additives : Passes as per FDA 21CFR 177.1520(a)1(i).

Particle Release : Passes test as per USP <788>, "Particulate matter in Injections".

Fiber Release : Complies with FDA 21CFR 210.3(b)(6).

CUSTOMER SUPPORT

mdi offers its unique interdisciplinary skills to provide solutions to specific problems. Please contact our factory or the local application specialist.

Chan

Head of Quality Assurance

Issue Date: 22-Jul-17

Advanced Microdevices Pvt. Ltd.

21, Industrial Area, Ambala Cantt, INDIA, Tel: +91-171-2699290/ 2699274 Website: www.mdimembrane.com Email: info@mdimembrane.com

An ISO 9001 Company