

# Certificate of Quality

The Microglassfiber syringe filters have been manufactured in a **mdi** facility in compliance with **ISO 9001** regulations.

## Microglassfiber Syringe Filters

Catalog No. : SYGF0605MNXX104  
Type : SY25GF  
Pore Size : 1.0 µm  
Lot Number : SG9587G

### SPECIFICATION

Membrane : Microglassfiber  
Housing : Polypropylene  
Filter Diameter : 25 mm  
Filtration Area : 4.15 cm<sup>2</sup>  
Burst pressure : > 14 Kg/cm<sup>2</sup>

### QUALITY ASSURANCE LOT RELEASE CRITERIA

Each lot has been sampled, tested and released by the Quality Assurance Department for the following characteristics:

- Integrity
- Flow Rate Performance
- Burst Pressure
- Leakage Test
- HPLC Certified

### BIOSAFETY

: Passes biological test for plastics Class VI as per USP and complies with Code of Federal Regulations FDA 21CFR 177.1520(a)1(i).

### EXPIRY

: Non degradable polymeric materials of construction with no expiry when stored at room temperature in a dry place.

### PRECAUTIONS

1. During handling, avoid contamination of outlet.
2. If pressure required to maintain the required flow becomes too high, the filter unit should be changed.

### CUSTOMER SUPPORT

**mdi** offers its unique interdisciplinary skills to provide solutions to specific problems. Please contact our factory or the local application specialist.



Head of Quality Assurance

Issue Date: 24-Jul-17

### Advanced Microdevices Pvt. Ltd.

21, Industrial Area, Ambala Cantt, INDIA,

Tel: +91-171-2699290/ 2699274

Website: [www.mdimembrane.com](http://www.mdimembrane.com)

Email: [info@mdimembrane.com](mailto:info@mdimembrane.com)

**An ISO 9001 Company**