

mdi Liner Drums with Trolleys

MDI cylindrical drum containers are robust, cost effective systems, specially designed and qualified for single use MDI *AseptiLiners* used for preparation, storage and transfer of media, buffer etc. in biopharmaceutical manufacturing. These plastic drums, constructed of high density polyethylene material, have smooth surfaces which are clean room compatible and are available with volumes of 50L, 100L, 200L and 500L.

The drums can be placed on specially designed MDI SS trolleys for safe and easy transfer.

MDI Liner Drums have applicability in large volume liquid handling operations such as storage, mixing and transfer.

Technical Specifications

Materials of construction

Drums : HDPE

Trolleys : Stainless Steel

Key Features

- > Ergonomic design for ease of use
- High Density Polyethylene drum construction
- Multi use application for shipping and storage
- Wheeled trolley available for portability
- 50L to 500L size option availability



Ordering Information

For Cylindrical Drums

Size	Catalog Number
50 L	PEDXHXXXXXXXX01
100 L	PEDXIXXXXXXXX01
200 L	PEDXJXXXXXXXX01
500 L	PEDXKXXXXXXXX01

For Wheeled SS Trolleys

Size	Catalog Number
50 L	TSSXHXXXXXXXX01
100 L	TSSXIXXXXXXXX01
200 L	TSSXJXXXXXXXX01
500 L	TSSXKXXXXXXXX01