

mdi BioKart SS trolley for 2D Bags

mdi BioKart is a specially designed stainless steel trolley with adjustable holder stand for 2D single use bags. These can hold 5L - 20L *AseptiBag* Gold 2D bag assemblies and are fully adjustable to a maximum height of 72 inches. The trolley has separate capsule filter holder to hold MDI large capsule filters of different sizes.

These trolleys support storage bag assemblies placed on the head and hang down from a stainless steel hook. The bottom of the trolley incorporates an twin compartment acrylic tray for holding accessories, bottles and tubing etc. The trolley is equipped with two front locking wheels for safety, when not in motion.

BioKart trolleys are made from high quality stainless steel construction which are easy to clean and disinfect. The mirror polished contact surface is smooth to ensure safety of the bags.

Key Features

- Ergonomic access
- > Corrosion resistant stainless steel construction
- > Enhanced mobility with strong wheels
- ➤ Lockable wheels for safety
- Convenient hanging
- Adjustable height to hand up to 20 liter bags
- > Twin compartment for segregated storage

Technical Specifications

Materials of construction

Trolley: SS316L **Tray**: Acrylic

Dimensions Height: 179 cm Length: 62.5 cm Width: 46 cm





Wheels with locking system



Clamp for height adjustment

Ordering Information

To order, please specify: BKSSTXXXXXXXXX1