

Certificate of Quality

The **ASESS** Aseptic Sampling Systems have been manufactured in a cGMP complaint **mdi** facility.

These products are produced using validated production processes in classified ISO 7 (Class 10, 000) area.

The quality management system complies with ISO 9001 standards.

ASESS Aseptic Sampling System

Catalog No. : A028XXXXXXSS301

Lot Number : AS0118E SI.No. 009

Ster. No. : R021

Date of Sterilization : 2018 - 05 Expiry Date : 2020 - 05

LOT RELEASE CRITERIA

The above lot of ASESS Sampling System meet the following lot release criteria:

100% Integrity Testing by pressure leak test

: Passes

100% Component Verification for specifications, orientation

and placement

: Conforms

100% Visual Inspection : Passes the visual inspection criteria

Sterilization : By gamma irradiation at ≥ 25 kGy

VALIDATED FOR

Bioburden Level : The bioburden level of sampling system components have been tested as per

ISO 11737

Sterility : The gamma sterilization process has been validated as per ISO 11137 to ensure

a Sterility assurance level (SAL) of 10⁻⁶

Endotoxin Level : < 0.25 EU/ml as determined by Limulus Amebocyte Lysate (LAL) test

Biosafety : Passes Biological Reactivity tests , In Vivo for Class VI plastic as per USP<88>

and Biological Reactivity tests, In Vitro for Cytotoxicity as per USP<87>

Particle release : Passes test as per USP<788> Particulate matter in Injections

Functionality: The sampling port of the sampling system withstands inline Steam Sterilization at

125 °C for 60 minutes and provides a sterile pathway for aseptic sampling

Shelf Life : Gamma irradiated ASESS Aseptic sampling system have been validated to have a

shelf life of 2 years

CUSTOMER SUPPORT

mdi offers its unique interdisciplinary skills to provide solutions to specific problems. Please contact our factory or the local application specialist.

Ouglity Assurance

Quality Assurance

Issue Date: 23-May-18

Advanced Microdevices Pvt. Ltd.

21, Industrial Area, Ambala Cantt, INDIA, Tel: +91-171-2699290/ 2699274 Website: www.mdimembrane.com Email: info@mdimembrane.com

An ISO 9001 Company