

Certificate of Quality

The Nylon-66 membrane syringe filters have been manufactured in a **mdi** facility in compliance with **ISO 9001** regulations.

Nylon-66 Membrane Syringe Filters

Catalog No. : SYNN0601MNXX204

Type : SY25NN-S
Pore Size : 0.2 μm
Lot Number : SN8898B

Ster. No. : 054
Expiry Date : 2021 - 02

SPECIFICATION

Membrane : Nylon-66
Housing : Polypropylene
Filter Diameter : 25 mm
Filtration Area : 4.15 cm²
Burst pressure :> 14 Kg/cm²
Hold-up volume :< 50 µl

Sterilization : Pre sterilized with EO

QUALITY ASSURANCE LOT RELEASE CRITERIA

Each lot has been sampled, tested and released by the Quality Assurance Department for the following characteristics:

- o Integrity by Bubble point Test using DI water Bubble point was: ≥ 50 psi (3.44 Bar)
- Flow Rate Performance
- o Burst Pressure
- o Leakage Test
- Retention: Microbial Challenge test as per ASTM F838-05 methodology

VALIDATED FOR

o Sterility

Blosafety : Passes Biological Reactivity Tests, In Vivo for Class VI plastic as described

in USP <88>.

Indirect Food Additives : Passes as per FDA 21CFR 177.1520(a)1(i).

PRECAUTIONS

- 1. During handling, avoid contamination of outlet.
- 2. If pressure required to maintain the required flow becomes too high, the filter unit should be changed.

CUSTOMER SUPPORT

mdi offers its unique interdisciplinary skills to provide solutions to specific problems. Please contact our factory or the local application specialist.

Our

Head of Quality Assurance

Issue Date: 20-Feb-18

Advanced Microdevices Pvt. Ltd.

21, Industrial Area, Ambala Cantt, INDIA, Tel: +91-171-2699290/ 2699274 Website: www.mdimembrane.com Email: info@mdimembrane.com

An ISO 9001 Company