

Certificate of Quality

The Nylon-66 membrane syringe filters have been manufactured in a **mdi** facility in compliance with **ISO 9001** regulations.

Nylon-66 Membrane Syringe Filters

| Catalog No. | : SYGN0601MNXX104 |
|-------------|-------------------|
| Туре | : SY25GN |
| Pore Size | : 0.2 µm |
| Lot Number | : SN7797G |

SPECIFICATION

| Membrane | : Nylon-66 with Multilayer Glassfiber Prefilter |
|-----------------|---|
| Housing | : Polypropylene |
| Filter Diameter | : 25 mm |
| Filtration Area | : 4.15 cm ² |
| Burst pressure | : > 14 Kg/cm² |

QUALITY ASSURANCE LOT RELEASE CRITERIA

Each lot has been sampled, tested and released by the Quality Assurance Department for the following characteristics:

- Integrity by Bubble point Test using DI water Bubble point was: ≥ 50 psi (3.44 Bar)
- Flow Rate Performance
- o Burst Pressure
- Leakage Test
- HPLC Testing

Biosafety : Passes Biological Reactivity Tests, *In Vivo* for Class VI plastic as described in USP <88>.

Indirect Food Additives : Passes as per FDA 21CFR 177.1520(a)1(i).

EXPIRY : Non degradable polymeric materials of construction with no expiry when stored at room temperature in a dry place.

PRECAUTIONS

- 1. During handling, avoid contamination of outlet.
- 2. If pressure required to maintain the required flow becomes too high, the filter unit should be changed.

CUSTOMER SUPPORT

mdi offers its unique interdisciplinary skills to provide solutions to specific problems. Please contact our factory or the local application specialist.

N-

Head of Quality Assurance Issue Date: 24-Jul-17

Advanced Microdevices Pvt. Ltd.

21, Industrial Area, Ambala Cantt, INDIA, Tel: +91-171-2699290/2699274 Website: www.mdimembrane.com Email: info@mdimembrane.com

An ISO 9001 Company