

Certificate of Quality

The Hydrophilic Polypropylene syringe filters have been manufactured in a **mdi** facility in compliance with **ISO 9001** regulations.

Hydrophilic Polypropylene Syringe Filters

Catalog No. : SYGP0602MNXX104

Type : SY25GP Pore Size : $0.45 \mu m$ Lot Number : SP9917G

SPECIFICATION

Membrane : Hydrophilic Polypropylene with Multilayer Glassfiber Prefilter

Housing : Polypropylene

Filter Diameter : 25 mm

Filtration Area : 4.15 cm^2 Burst pressure : > 14 Kg/cm^2 Hold up volume : < 50 µl

QUALITY ASSURANCE LOT RELEASE CRITERIA

Each lot has been sampled, tested and released by the Quality Assurance Department for the following characteristics:

- o Integrity
- o Flow Rate Performance
- Burst PressureLeakage Test
- HPLC Test
- **BIOSAFETY**: Passes biological test for plastics Class VI as per USP and complies with Code of

Federal Regulations FDA 21CFR 177.1520(a)1(i).

EXPIRY : Non degradable polymeric materials of construction with no expiry when stored at

room temperature in a dry place.

PRECAUTIONS

- 1. During handling, avoid contamination of outlet.
- 2. If pressure required to maintain the required flow becomes too high, the filter unit should be changed.

CUSTOMER SUPPORT

mdi offers its unique interdisciplinary skills to provide solutions to specific problems. Please contact our factory or the local application specialist.

Pin-

Head of Quality Assurance

Issue Date: 24-Jul-17

Advanced Microdevices Pvt. Ltd.

21, Industrial Area, Ambala Cantt, INDIA, Tel: +91-171-2699290/ 2699274 Website: www.mdimembrane.com Email: info@mdimembrane.com

An ISO 9001 Company