

Certificate of Quality

The Nylon-66 membrane syringe filters have been manufactured in a **mdi** facility in compliance with **ISO 9001** regulations.

Nylon-66 Membrane Syringe Filters

Catalog No. : SYNN0602MNXX102

Type : SY25NN
Pore Size : $0.45 \mu m$
Lot Number : SN5567G

SPECIFICATION

Membrane : Nylon-66
Housing : Polypropylene
Filter Diameter : 25 mm
Filtration Area : 4.15 cm²
Burst pressure :> 14 Kg/cm²
Hold-up volume :< 50 µl

QUALITY ASSURANCE LOT RELEASE CRITERIA

Each lot has been sampled, tested and released by the Quality Assurance Department for the following characteristics:

Integrity

Flow Rate Performance

Burst PressureLeakage TestHPLC Testing

Biosafety : Passes Biological Reactivity Tests, In Vivo for Class VI plastic as described

in USP <88>.

Indirect Food Additives : Passes as per FDA 21CFR 177.1520(a)1(i).

EXPIRY : Non degradable polymeric materials of construction with no expiry when stored at

room temperature in a dry place.

PRECAUTIONS

1. During handling, avoid contamination of outlet.

2. If pressure required to maintain the required flow becomes too high, the filter unit should be changed.

CUSTOMER SUPPORT

mdi offers its unique interdisciplinary skills to provide solutions to specific problems. Please contact our factory or the local application specialist.

Com

Head of Quality Assurance

Issue Date: 24-Jul-17

Advanced Microdevices Pvt. Ltd.

21, Industrial Area, Ambala Cantt, INDIA, Tel: +91-171-2699290/ 2699274
Website: www.mdimembrane.com
Email: info@mdimembrane.com

An ISO 9001 Company