

## PTFE Membrane Syringe Filter with Microglassfiber Pre-filter Type- SYTG

**mdi** SYTG syringe filters with PTFE membrane and multilayered binder free borosilicate microglassfiber pre-filter provide absolute retention of very difficult to filter colloidal contaminants.

These are specially used to filter highly turbid (hazy) solutions.

## **Special Features**

- Absolute retention efficiency
- 16 Channels for minimum back pressure
- Extremely low extractable
- Wide chemical compatibility
- HPLC certified for low levels of UV absorbing extrac-tables
- Heat sealed to ensure 'no leaching'
- No initial flushing required

#### **Applications**

- Sample filtration for HPLC and other critical instruments
- Filtration of organic solvents and chemicals





## Specifications

#### Materials of Construction

Housing	Polypropylene
Pre-filter	Microglassfiber
Filter Media	PTFE Membrane

**Pore Size** 0.2 μm, 0.45 μm

#### Inlet/Outlet Connections Inlet: Female Iuer-lock Outlet: Male Iuer slip

# Typical Hold-up Volume(with air purge)

13mm	25mm
<20µl	<50µl

Burst Pressure (13mm/25mm) >14kg/cm<sup>2</sup>

Maximum Operating Temperature  $55\ ^{\circ}\mathrm{C}$ 

Maximum Operating Pressure 75 psi @ 25° C

**Recommended Sample Volume** 

13mm 25mm

< 10 ml < 100 ml



## Fiber Release

Complies with USFDA CFR Title 21, Part 210.3(b)(6) and do not requires any initial Flushing.

#### **Particle Release**

The filtrate Complies with USP <788> test for particulate matter in injections.

#### **Indirect Food Additives**

The plastic housing material complies with the requirements of USFDA CFR Title 21 Part 177.1520

#### **HPLC Certified**

SYTG syringe filters are tested and certified for use in HPLC applications

#### **Oxidizable Matter**

Passes test as per USP <1231>

DST SYTGXXX1248F

Data Sheel



## Packaging

### Standard Pack:

SYTG syringe filters are packed in PET jars, duly labeled. Each jar has 100 Nos. syringe filters.

#### **Bulk Pack:**

SYTG syringe filters are also bulkpacked in polyethene bags, duly labeled. Each bag has 1000 Nos. syringe filters.

## Storage

i.

It is recommended to keep the syringe filters:

- at clean and dry location
- ii. away from direct heat sources at temperature below 50°C

## **Ordering Information**

Туре		Size		Pore Size		Inlet/Outlet		ХХ	Sterility		Pack Size	
	Code	Diameter	Code		Code		Code			Code	Qty	Code
PTFE Membrane Syringe Filter with SY microglassfiber pre-filter	SYTG	13mm	03	0.2µm	01	Female Luer Lock	М		Non Sterile	1	100	04
	5116	25mm	06	0.45µm	02	Male Luer Slip	N				1000	06

SYTG 06 01 MN XX 1 04	Exa	imple:						
		SYTG	06	01	MN	XX	1	04