

Certificate of Quality

The PVDF membrane syringe filters have been manufactured in a mdi facility in compliance with ISO 9001 regulations.

PVDF Membrane Syringe Filters

: SYVF0601MNXX104 Catalog No.

Type : SY25VF Pore Size : 0.2 µm Lot Number : SV8467G

SPECIFICATION

: PVDF Membrane

Housing : Polypropylene Filter Diameter : 25 mm Filtration Area : 4.15 cm² Burst pressure : > 14 Kg/cm²

QUALITY ASSURANCE LOT RELEASE CRITERIA

Each lot has been sampled, tested and released by the Quality Assurance Department for the following characteristics:

- Integrity by Bubble point Test using DI water Bubble point was: ≥ 45 psi (3.10 Bar)
- Flow Rate Performance
- **Burst Pressure** Leakage Test
- Retention: Microbial Challenge test as per ASTM F838-05 methodology
- **HPLC** Testing

Biosafety : Passes Biological Reactivity Tests, In Vivo for Class VI plastic as described

in USP <88>.

Indirect Food Additives : Passes as per FDA 21CFR 177.1520(a)1(i).

EXPIRY : Non degradable polymeric materials of construction with no expiry when stored at

room temperature in a dry place.

PRECAUTIONS

1. During handling, avoid contamination of outlet.

2. If pressure required to maintain the required flow becomes too high, the filter unit should be changed.

CUSTOMER SUPPORT

mdi offers its unique interdisciplinary skills to provide solutions to specific problems. Please contact our factory or the local application specialist.

21, Industrial Area, Ambala Cantt, INDIA, +91-171-2699290/ 2699274 Website: www.mdimembrane.com Head of Quality Assurance info@mdimembrane.com

An ISO 9001 Company

Advanced Microdevices Pvt. Ltd.

Issue Date: 24-Jul-17