

Microglassfiber Prefilter Type- GF2 (For High Contamination Load)

mdi Microglassfiber Prefilter type-GF2 are high dirt holding capacity filters most widely used as a pre-filter for membranes and to clarify various solutions in the laboratory.

mdi filters are produced in a facility with ISO-9001:2008 certified quality management system.



Special Features

- Good retention efficiency
- High flow rates
- High dirt holding capacity
- Wide chemical compatibility
- Biologically inert

Applications

- Pre-filtration of solutions with high dirt load

Specifications

Materials of Construction

Filter	Borosilicate
Media	Microglassfiber

Size

10mm, 13mm, 25mm, 47mm, 90mm, 127mm, 142mm, 257mm, 279mm, 293mm

Pore Size

1.5µm

Thickness

450 ± 50µm

Extractables with WFI

Passes test as per USP

Oxidizable Substances

Within limits as specified in USP

Ordering Information

Type		Size		Pore Size		XX	XX	Sterility		Pack Size	
	Code	Diameter	Code		Code				Code	Qty	Code
GF2	GF2X	10mm	02	1.5µm	14			Non Sterile	1	100	04
		13mm	03								
		25mm	06								
		47mm	09								
		90mm	14								
		127mm	15								
		142mm	16								
		257mm	17								
		279mm	18								
		293mm	19								
Example:											
GF2X		09		14		XX	XX	1			04