

Certificate of Quality

The Nylon-66 membrane centrifugal filters have been manufactured in a **mdi** facility in compliance with **ISO 9001** regulations.

Pre-sterilized Nylon-66 Membrane Centrifugal Filters

Catalog No. : CFNN2101XXXX203

Type : CF7NN-S

Pore Size : 0.2 μm

Lot Number : FN8877G

Ster. No. : 025

Expiry Date : 2020 - 07

SPECIFICATION

Membrane : Nylon-66 Housing : Polypropylene

Filter Diameter : 7mm
Filtration Area : 0.28cm²
Hold-up volume : < 5µl
Operating pH range : 4 to 14
Maximum sample volume : 600µl

Maximum operating pressure : 5600 x g at 10,000 rpm

Maximum operating temperature : 80 °C max

Sterilization : Pre sterilized with EO

100% Integrity Tested : Each unit is tested for integrity during manufacture.

QUALITY ASSURANCE LOT RELEASE CRITERIA

Each lot has been sampled, tested and released by the Quality Assurance Department for the following characteristics:

- o Integrity
- Flow Rate Performance
- HPLC Certified
- Sterility

CUSTOMER SUPPORT

mdi offers its unique interdisciplinary skills to provide solutions to specific problems. Please contact our factory or the local application specialist.

Com

Head of Quality Assurance

Issue Date: 24-Jul-17

Advanced Microdevices Pvt. Ltd.

21, Industrial Area, Ambala Cantt, INDIA, Tel: +91-171-2699290/ 2699274
Website: www.mdimembrane.com
Email: info@mdimembrane.com

An ISO 9001 Company