# Filters For Microbiology



# Bio-burden testing of process water and finished products in

- Pharmaceuticals
- Biopharmaceuticals
- Food & Beverages
- Cosmetics







# M-Funnels

## **EO Sterilized, Ready to Use**



#### **Variants**

- With Media Reservoir Pad
- Without Pad

#### Sizes Available

➤ 100ml
➤ 250ml

#### Validated for

- Microbial Recovery
  - With Liquid Media
  - With Semi-solid Agar Media
- > ASTM 4200-82 (2003)
  - Evaluating inhibitory effects of ink grids on membrane filters
- Microbial Retention
- Sterility



# Unique Performance Advantage

#### **Engraved volume markings**

For accurate measurement of sample volumes



#### Validated ink Grids

- No inhibitory effect on microbial recovery (ASTM 4200-82)
- No inference with colony counting



# **Universally Adaptable Design**

Fits directly on most filter manifolds

### Validated "Media Reservoir" Pad\*

For use with 2ml liquid media





#### Special Base Design

For safe and easy lifting of membrane

#### Easy to Incubate Design\*

Converts into a petri plate with an easy to "Squeeze Remove" cup and transparent lid



# **Ordering Information**

Туре		Si	ze	Pore S	ize	Х	Support	Pad	Capa	city	Sterilizat	ion	Pack S	Size
Type	Code	Dia	Code	Pore Size	Code			Code		Code		Code	Pack Size	Code
M-Funnel	FMCN	47mm	09	0.45µm	02		without pad	Х	100 ml	XX	EO Sterile	2	24 (100 ml)	12
•							with pad	Р	250 ml	01			12 (250 ml)	08

#### **Example:**

FMCN	09	02	X	Р	01	2	08
	1			I			

<sup>\*</sup>Applicable for 'with pad' variant only

# Grid Marked Cellulose Nitrate Membrane Disc Filters

## **Type Available**

- > GCN
- > RGCN

#### Validated for

- Microbial Recovery
- > ASTM 4200-82 (2003)
  - Evaluating inhibitory effects of ink grids on membrane filters
- Microbial Retention
- Sterility



## **Ordering Information**

Туре						
Type	Code					
GCN	GCNX					
RGCN RGCNX						

Size					
Dia	Code				
47mm	09				
50mm	10				

Pore Size							
Pore Size	Code						
0.45µm	02						

хх	

хх	Sterile /Non Sterile				
		Cod			
	Non Sterile	1			
	EO Sterile	2			

Sterile Non Ste	-		Pack S	Size
Code			Pack Size	Code
on Sterile	1		100	04
O Sterile	2		150	15

Note:

GCN is grid marked Cellulose Nitrate membrane available in non-sterile as well as individually packed pre-sterilized packs. Available in pack size of 100 Nos only

RGCN is individually packed pre sterilized GCN membrane in reel form. Available in pack size of 150 Nos only

**Example:** RGCN 09 02 XX XX 2 15

#### mdi Membrane Technologies INC.

5340 Jaycee Avenue, Suite A, Harrisburg, PA 17112 Phone: +1-717-412-0943, Fax: +1-717-695-9637 E-mail: rs@mdimembranetech.com Website: www.mdimembranetech.com 20-21, Industrial Area, Ambala Cantt- 133 006, India Tel: +91 - 171-2699290, 2699471 E-mail: info@mdimembrane.com Website: www.mdimembrane.com