

Certificate of Quality

The Polyethersulfone Membrane Vacuum Filter Device has been manufactured in a **mdi** facility in compliance with **ISO 9001** regulations using *validated production processes*.

AseptiVac KL-y PES Membrane Vacuum Filter Device

Catalog No. : AKLX3101D10X302

Type : AVKL-S Pore Size : 0.2 μm

Lot Number : AK8827G

Ster. No. : R011

Expiry Date : 2019 - 07

SPECIFICATIONS

Materials of Construction

Filter Membrane : Polyethersulfone
Filter Housing : Polypropylene
Receiver : Polystyrene
Filtration Area (Nominal) : 100 cm²
Maximum Operating Temperature : 45 °C

Sterilization : Pre sterilized by Gamma Irradiation

LOT RELEASE CRITERIA

100% Integrity Tested : The filter has been tested for integrity by Bubble Point Test using 50% IPA/Water

solution. Bubble point was: ≥ 16 psi (1.10 Bar)

Water Flow Rate : ≥ 1.4 lpm @ Vacuum of 500 mm of Hg @ 27 °C

Microbial Challenge Test : Retains $\geq 10^7$ organisms/cm² of *B. diminuta* ATCC 19146 challenge as per ASTM

F838-05 methodology.

Sterility : Samples passed the sterility test in accordance with USP <71>.

VALIDATED FOR

Bacterial Endotoxin Levels : Aqueous extracts exhibit < 0.25 EU/ml as established by Limulus Amebocyte Lysate

(LAL) test.

Biosafety : Passes Biological Reactivity Tests, In Vivo for Class VI plastic as described in USP

<88>.

Cytotoxicity: Passes Biological Reactivity Tests, *In Vitro* as described in USP <87>.

Extractable: Within limits as specified in USP.

CUSTOMER SUPPORT

mdi offers its unique interdisciplinary skills to provide solutions to specific problems. Please contact our factory or the local application specialist.

Que

Head of Quality Assurance

Issue Date: 24-Jul-17

Advanced Microdevices Pvt. Ltd.

21, Industrial Area, Ambala Cantt, INDIA, Tel: +91-171-2699290/ 2699274 Website: www.mdimembrane.com Email: info@mdimembrane.com

An ISO 9001 Company