

Certificate of Quality

The Absorbent Pad has been manufactured in a **mdi** facility in compliance with ISO-9001 regulations using validated production processes.

ABSORBENT PAD

TYPE : AP120

LOT NO. : S161-5274

PHYSICAL TESTS	VALUES	SPECIFICATION
Thickness (μm)	: 1491	1350 – 1800 μm
Tare weight (mg/cm^2)	: 78.4	70 – 95 mg/cm^2
Width/Dimensions (mm)	: 20x7.5	20 x 7.5 (± 0.5) mm

The dimensions are measured after every 5000 cuts and have been found in specification : Length - 20 ± 0.5 mm, Height – 7.5 ± 0.5 mm and Angle $80^\circ \pm 3^\circ$.



Head of Quality Control

Issue Date: 04-Jul-2013

Advanced Microdevices Pvt. L
21, Industrial Area, Ambala Cantt, IN
Tel: +91-171-2699290/2699
Fax: +91-171-2699221/2699
Website: www.mdimembrane.c
Email: mdi@vsnl.c

An ISO 9001 Compa