

Rapid Test Components Cover Tapes

Product CoverTapes

Intended Use Manufacture of Lateral Flow Immunoassay Test Kit. The Clear Tape is used to cover

the membrane so that the membrane does not dry up and at the same time the user can see the lines that develop. The Masking Taps are used to hold the pads and

membrane together with the plastic backing material.

Type Masking Tape (MT1 and MT2)

 $Clear Tape\,CT$

Physical Properties

Property	MT1	MT2	СТ	
Appearance	White	White	Transparent	
Thickness	180 - 280 μm		80 - 140 μm	
Adhesion	The tape does not comes off the membrane plastic after immersion in water for 24 hours.			

Supply Format

	Masking Tape	Clear Tape	
Format	Strip	Strip/Reel	
Dimensions	Desired Dimensions		

Storage Condition Store and handle in accordance with all current regulations and standards.

Store in original packing at 15-30 °C. Avoid direct sunlight.

Caution: Material should be kept away from heat radiators.

Shelf Life One Year from the date of dispatch, stored under above conditions.

Ordering Information						
Product		Description				
Туре		Strip/Reels Strip/Reels Length Width		Quantity (Nos.)		
*MT1		As per requirement	**As per requirement	As per requirement		
*MT2						
СТ		* MT1 and MT2 is only available in strip format ** Maximum strip/reels width is 5cm				
xample: ♥	\downarrow	\downarrow	↓	\downarrow		
MT2		30 cm	3 cm	-		