

# Certificate of Quality

The Conjugate Release Matrix has been manufactured in a **mdi** facility in compliance with ISO-9001 regulations using validated production processes.

## CONJUGATE RELEASE MATRIX

**TYPE** : PTR-5

**LOT NO.** : 681813-75

## PHYSICAL TESTS

Thickness ( $\mu\text{m}$ ) : 404  
Tare weight ( $\text{mg}/\text{cm}^2$ ) : 9.8  
H<sub>2</sub>O Wicking Rate, 4cm distance (sec) : 24  
Width (mm) : 8

## FUNCTIONAL TESTS

Accelerated Gold Release Test : No visible mark  
45 °C/24 hours : at point of application

**BEST BEFORE** : One year from the date of dispatch.

## STORAGE

For best results Conjugate Release Matrix should be stored in original packing at 15 – 30 °C temperature. It should be kept away from direct sun light and heat radiators.



Head of Quality Control

Issue Date: 04-Jul-2013

**Advanced Microdevices Pvt. Ltd.**

21, Industrial Area, Ambala Cantt, INDIA,

Tel: +91-171-2699290/2699274

Fax: +91-171-2699221/2699008

Website: [www.mdimembrane.com](http://www.mdimembrane.com)

Email: [mdi@vsnl.com](mailto:mdi@vsnl.com)

**An ISO 9001 Company**