

# **Certificate of Quality**

The pDNA Miniprep Kits have been manufactured in mdi facility in compliance with ISO 9001 regulations using validated production processes.

## **pDNA Miniprep Kit**

Catalog No. : MIPKXXXXXXX0250

Type : MIPK

Lot Number : PD1-1693-122

#### **LOT RELEASE CRITERIA**

Binding capacity of column : > 25  $\mu$ g pDNA Yield \* : up to 25  $\mu$ g DNA Purity (A<sub>260</sub> /A<sub>280</sub>) : 1.8 – 2.0 DNase Free : Yes

Plasmid type – pUC vectors \*
Copy number - High (500-700)

Fragment size - 9 kb Culture volume – up to 5ml Elution volume – 50µl

#### **VALIDATED FOR**

- Downstream Processing
  - Restriction Digestion
  - Ligation
  - PCR
  - Sequencing
  - Bacterial Transformation
  - Cloning
- Shelf Life
  - Minimum one year at Room Temperature
  - More than two years at 2-8°C
- Purification of low copy number plasmids as well as high copy number plasmids (<10 copies per cell)
- Purification of small plasmids (2 kb) as well as large plasmids (10 kb)
- Consistent intra lot (column to column) and inter lot plasmid DNA recovery

### **CUSTOMER SUPPORT**

**mdi** offers its unique interdisciplinary skills to provide solutions to specific problems. Please contact our factory or the local application specialist.

Head of Quality Control Issue Date: 28-Jun-13

Tel: +91-171-2699290/ 2699274
Fax: +91-171-2699221/ 2699008
Website: www.mdimembrane.com
Email: mdi @vsnl.com

**Advanced Microdevices Pvt. Ltd.** 21, Industrial Area, Ambala Cantt, INDIA,

. ... ....

An ISO 9001 Company