

Certificate of Quality

The **Endotoxin free Quanta Mega Kits** have been manufactured in **mdi** facility in compliance with **ISO 9001** regulations using **validated production processes**.

Endotoxin free Quanta Mega Kit

Catalog No. : QMEKXXXXXXXXX0002
Type : QMEK
Lot Number : QM1-1473-107

LOT RELEASE CRITERIA

Binding capacity of Column : ≥ 3 mg
***pDNA Yield** : Upto 2.5 mg
DNA Purity (A_{260}/A_{280}) : 1.8 – 2.0
DNase Free : Yes
Endotoxin Free : ≤ 0.1 EU/ μ g

*Plasmid type – pUC

Copy number – High (500-700)

Fragment Size – 9 kb

Culture Volume– 2.5 liters

Elution volume – 1 ml

VALIDATED FOR

- Downstream Processing
 - Restriction Digestion
 - Transfection
 - PCR
 - Sequencing
 - Bacterial Transformation
 - Cloning
- Shelf Life
 - Minimum two years at Room Temperature
- Purification of low copy number plasmids (< 10 copies per cell) as well as high copy number plasmids.
- Purification of small plasmids as (2 kb) as well as large plasmids.
- Consistent intra lot (column to column) and inter lot plasmid DNA recovery

CUSTOMER SUPPORT

mdi offers its unique interdisciplinary skills to provide solutions to specific problems. Please contact our factory or the local application specialist.



Head of Quality Control

Issue Date: 28-Jun-13

Advanced Microdevices Pvt. Ltd.

21, Industrial Area, Ambala Cantt, INDIA,

Tel: +91-171-2699290/ 2699274

Fax: +91-171-2699221/ 2699008

Website: www.mdimembrane.com

Email: mdi@vsnl.com

An ISO 9001 Company