

Certificate of Quality

The Endotoxin Free Quanta Midi Kits have been manufactured in mdi facility in compliance with ISO 9001 regulations using validated production processes.

Endotoxin Free Quanta Midi Kit

Catalog No. : QDEKXXXXXXX0002

Type : QDEK

Lot Number : QD1-1453-108

LOT RELEASE CRITERIA

Binding capacity of Column : $\geq 350 \ \mu g$ *pDNA Yield : Up to 250 μg DNA Purity (A₂₆₀ /A₂₈₀) : 1.8 – 2.0 DNase Free : Yes

Endotoxin Free : ≤ 0.1 EU/µg

*Plasmid type - pUC

Copy number – High (500-700)

Fragment Size – 9 kb Culture Volume– 35 ml Elution volume – 200 µl

VALIDATED FOR

- · Downstream Processing
 - > Restriction Digestion
 - > Transfection
 - ➤ PCR
 - Sequencing
 - Bacterial Transformation
 - Cloning
- Shelf Life
 - Minimum two years at Room Temperature
- Purification of low copy number plasmids as well as high copy number plasmids (< 10 copies per cell)
- Purification of small plasmids as well as large plasmids (2 kb)
- · Consistent intra lot (column to column) and inter lot plasmid DNA recovery

CUSTOMER SUPPORT

mdi offers its unique interdisciplinary skills to provide solutions to specific problems. Please contact our factory or the local application specialist.

Head of Quality Control Issue Date: 28-Jun-13

Advanced Microdevices Pvt. Ltd.

21, Industrial Area, Ambala Cantt, INDIA, Tel: +91-171-2699290/ 2699274 Fax: +91-171-2699221/ 2699008 Website: www.mdimembrane.com

Email: mdi @vsnl.com

An ISO 9001 Company