



Membrane Technologies

Steridisc Dispenser Model SD-01

mdj Steridisc Dispenser is an advanced, light weight table top instrument for dispensing individually packed membrane disc filters from the pre-sterilized reels.

This all-new fully automatic dispenser with electronic control is very easy to operate. The push-button system makes disc handling very convenient for the user. It delivers a single sterile membrane and simultaneously wraps the peeled-off protective packaging.

The stainless steel clips hold the free end of the packaging material (paper on one side and the plastic film on the other side), while releasing the membrane disc filter, absolutely untouched by hand.

A special feature of the Steridisc Dispenser is operation by foot paddle, which allows totally hands free operation for the microbiologist.

Specifications

Dimensions

Height	: 270 mm
Width	: 175 mm
Length	: 250 mm

Weight : 3 Kg

Voltage : 100-240V AC

Frequency : 50 Hz

Dispense Speed : < 2 seconds per disc (fixed)

Storage : It is recommended to be stored in a clean environment at 10° to 40° C

Applications

- Membrane dispensing for sterility testing and microbiology



Advantages

- Automatic single unit dispensing
- Less handling time
- No accidental contamination
- Easy loading of membrane roll
- Easy portability: Small and light
- Push-button or foot paddle operated dispensing
- Easy to clean and maintain

Ordering Information

Type	
mdj Steridisc Dispenser	Model SD-01

To Order please Specify:		
Type	Model No.	Qty.