

## mdi BioSafe

# Stainless Steel Rectangular Totes

**mdi** BioSafe stainless steel rectangular totes are robust systems, specially designed for secure transportation of **mdi** AseptiBag Gold 3D bags filled with buffers, media, process intermediates and drug substances in any stage of bio manufacturing process within or outside the facility.

These rectangular tote boxes are made from high quality stainless steel construction which are easy to clean and disinfect and the smooth mirror polished contact surface ensures bag safety.

Stainless steel totes ranging from 100 liters to 500 liters have applicability in large volume liquid-handling operations such as storage, transfer and waste collection.

## **Application**

- Storage and transportation of buffers, media, process intermediates and drug substances
- > Waste collection

### **Key Features**

- > Ergonomic design
- > 304 L stainless steel corrosion resistant construction
- Smooth mirror polished contact surface to ensures bag safety
- > Enhanced mobility with strong wheels
- > Front brakes for security
- Large door for easy access and positioning of process containers
- > SS latch on door to prevent accidental opening
- > Bottom drain port for easy installation of drain tubing

## Ordering Information

With Wheels

For 100L/200 L bag : RTSSJXXXXWXXX01
For 500 L bag : RTSSKXXXXWXXX01



## **Technical Specifications**

#### **Materials of construction**

304 L Stainless Steel

#### **Surface Finish**

Mirror Polished

### **Dimensions and Weight**

Tote Box	Internal Dimensions (LXBXH) mm	External Dimensions (LXBXH) mm	Weight (kg)
100 L	730 x 530 x 575	810 x 610 x 1000	76
200 L	730 x 530 x 575	810 x 610 x 1000	76
500 L	1100 x 720 x 730	1180 x 800 x 1120	110

#### **Without Wheels**

For 100L/200 L bag : RTSSJXXXXXXXXX01
For 500 L bag : RTSSKXXXXXXXXX01