

Certificate of Quality

The Genomic DNA Miniprep Kits have been manufactured in **mdi** facility in compliance with ISO 9001 regulations using *validated production processes*.

Genomic DNA Miniprep Kit

Catalog No.

Туре

- Lot Number
- : CTGKXXXXXXX0050
- : CTGK
- : CT1-0660-040

LOT RELEASE CRITERIA

Binding capacity of column gDNA Yields *

:	>	50	μg
---	---	----	----

:	Sample Type	Yields	
	Mammalian tissue	Upto 30 µg	
	Bacterial cells	Upto 22 µg	
	Mammalian blood	Upto 12 µg	

DNA Purity (A₂₆₀ /A₂₈₀) DNase Free

: 1.7 – 1.9

: Yes

Amount of tissue-25 mg* Amount of Bacterial Cells-2 X 10⁹ Amount of Mammalian blood-200 μl Elution volume – 200 μl

VALIDATED FOR

- Downstream Processing
 - Restriction Digestion
 - PCR
 - Sequencing
 - Cloning
- Shelf Life
 - > Minimum one year at Room Temperature
 - More than two years at 2-8°C
- Purification of small as well as large gDNA fragments
- Consistent intra lot (column to column) and inter lot genomic DNA recovery

CUSTOMER SUPPORT

mdi offers its unique interdisciplinary skills to provide solutions to specific problems. Please contact our factory or the local application specialist.

Head of Quality Control

Issue Date: 23-Aug-11

Advanced Microdevices Pvt. Ltd.

 21, Industrial Area, Ambala Cantt, INDIA,

 Tel:
 +91-171-2699290/2699274

 Fax:
 +91-171-2699221/2699008

 Website:
 www.mdimembrane.com

 Email:
 mdi@vsnl.com

An ISO 9001:2000 Company